




you expect
Ease
you demand
Performance

QuickTake 30 Sample Pump


- **Grab sample pump with constant flow from 10 to 30 L/min**
 - Maintains set flow throughout a sampling period for complete sample integrity
 - Ideal for use with BioStage Impactors, spore trap cassettes such as VersaTrap or Air-O-Cell, asbestos cassettes, wall cavity samplers, or impaction samplers that require flows to 30 L/min
- **Run times up to four hours on one battery charge**
 - AC charger/adaptor for extended sampling
 - Reliable rechargeable lithium-ion battery pack
- **Compact and portable**
 - 9.3 x 8.4 x 3.5 in (23.6 x 21.3 x 8.9 cm)
 - Weighs only 4.8 pounds (2.2 kg)
 - Handle for easy portability
- **Easy-to-use programmable timer for unattended sampling**
 - Select sample time presets or customize timer from 1 to 999 minutes
 - Continuous run with manual shut-off
 - Intermittent sampling
 - Large, easy-to-use keypad
- **Versatile sampling accessories**
 - Wand for sampling in ductwork and other hard-to-reach areas
 - Mounting bracket for BioStage allows impactor to be mounted directly on pump
- **Easy-to-use**
 - Simple no-tool calibration using any reliable calibrator
 - End-of-cycle visual and audible alarms



QuickTake 30 with BioStage on mounting bracket



QuickTake 30 with spore trap cassette



QuickTake 30 with spore trap cassette on sampling wand

QuickTake 30 Sample Pump Ordering

QuickTake Pumps contain Li-Ion batteries and are subject to special shipping regulations.

Description	Cat. No.
QuickTake 30 Sample Pump includes battery pack, AC charger/adaptor (100 - 240 V), and cassette and tubing adapters	228-9530
QuickTake 30 BioStage Pump Kit includes pump as described above, Standard BioStage, mounting bracket with inlet adapters, rotameter, & deluxe carry case	228-9530K
QuickTake 30 Sample Pump with Rotameter includes battery pack, AC charger/adaptor (110-240V), and cassette and tubing adapters	228-9530A
Accessories	
Mounting Bracket for BioStage/BioCassette includes inlet adapters	228-9531
Sampling Wand , telescopes to 34 ins for ductwork and hard-to-reach areas. Includes 6-ft tubing.	228-9521
Rotameter for QuickTake 15 and 30	320-100

