

Metals Wipe Media LOQ



LOQ µg

ANALYTES

	Lead Wipes ICP-MS	Lead Wipes ICP	Gauze Wipes ICP	Gauze Wipes ICP-MS
Aluminum			25	25
Antimony	3	3	3	3
Arsenic	1	2	1	1
Barium	5.0	5.0	5.0	5.0
Beryllium	0.05	0.5	0.5	0.1
Boron			5.0	
Cadmium	0.05	0.5	0.5	0.5
Chromium	7.5	7.5	5	5
Cobalt	1.5	1.5	1.5	1.5
Copper	4	4	1	1
Iron	250	250	50	50
Lead	1.3	1.3	1.3	1.3
Magnesium	500	250	250	250
Manganese	15	15	2.5	
Molybdenum	0.5	0.5	1.0	1
Nickel	2	2	1	1
Selenium	2.5	2.5	2.5	2.5
Tellurium				0.5
Thallium	2.5	5	5	5
Titanium	7	7	2.5	2.5
Vanadium	2.5	2.5	1.5	1.5
Zinc	200	200	25	25

For any other analytes, please call client services at 1-888-432-5227.